
## GETINTHE COMFORT ZONE

with new 3D printed custom orthotics

Thousands of data points precisely capture the shape and movement of your feet.



Advanced technology enables your foot specialist to design the perfect orthotic.

| GO-4D REPORT - LEFT FOOT  |   |
|---|---|
| Diagnostic  | Customization   |
| Too much pressure in the heel -40.3% on on left foot                      | Heel spur pad<br>+30.2%   |
| Overpronate (dropped arch)<br>-15.8% on on left foot                      | Raise the arch in height,<br>width and length<br>H 80mm, L 195mm, W 109mm |
| Excessive pressure and time on the ball of the foot -8.6% on on left foot | Metatarsal lift to spread out<br>the pressure and timing<br>+24.9%        |

the most corrective and comfortable fit you will ever experience.





Ask your footcare specialist about Go 4-D orthotics.